

How to Verify the FP40's Compatability with DSL 4.1

10/31/01

FP40 Instructions:

- 1) Select the "SPL" approach:
 - a) From the Main Screen press F5 to enter the "PROBE" screen.
 - b) Press F2 to select "SPL" approach.
- 2) Select the Fitting Rule:
 - a) Press F4 to enter the "TARGET" screen.
 - b) Use F3 to select "DSL LIN" or "DSL WDRC".
 - c) Press F5 to make this the current rule.
- 3) Select the patients age range using F6.
- 4) Select "HTL" with F2 and enter the thresholds into the "HTL" column.
- 5) If custom UCL is desired:
 - a) Select "UCL" with F2.
 - b) Enter the UCL's into the "UCL" column.
- 6) If DSL WDRC rule is used, select the Compression Threshold:
 - a) Select "CT" with F2.
 - b) Dial in the Compression Threshold with the AMPLITUDE control.
- 7) Select the Aid Type:
 - a) Press F7 to enter the "COUPLER" screen.
 - b) Use F4 to select the aid type.
- 8) Select the RECD:
 - a) Press F3 to enter "RECD MODE".
 - b) If average RECD is desired, select "AVG" RECD with F3.
 - c) If custom RECD is desired:
 - i) Enter the RECD's into the RECD column.
 - ii) Select "MEASURED RECD" with F3.
- 9) Generate the Targets:
 - a) Press F7 to return to the "TARGET" screen.
 - b) Press F5 to generate the targets.

How to Verify the FP40's Compatability with DSL 4.1

10/31/01

FP40 Instructions(cont):

- 10) View the Real Ear Speech-Weighted Input targets:
 - a) Press F4 to enter the "PROBE" screen.
 - b) Make sure the source in the curve box for CRV 2 is either "70 COMP" or "60 NORM".
 - c) If necessary, measure and then undisplay curve 2 to make it so.
 - d) The targets for "T"reshold, "A"verage Speech, and "U"CL's can be viewed as character graphs.
 - e) The Average Speech target can also be viewed as a tabulation.
 - i) Press F1 to enter the Setup Menu.
 - ii) Select DATA DISPLAY - TARG SPL under PROBE SETTINGS.
 - iii) Press F1 to return to the Probe screen and press the DATA/GRAPH key.
- 11) View the Coupler Gain Average Speech Weighted Input Target:
 - a) Measure a Composite curve in the "PROBE" screen.
 - b) Press F4 and then F7 to enter the "COUPLER" screen.
 - c) Make sure the "MODS" and "RESER" are all 0's.
 - d) Use F6 to set the source to 70.
 - e) The Coupler Gain Target can now be viewed. Use DATA/GRAPH to select graph or tabular form.
 - f) The Compression Ratio's for each freq can also be viewed if the DSL WDRC rule is used.
- 12) View the Coupler Gain Average Constant Level Input Target:
 - a) Press F7 and then F4 to return to the "PROBE" screen.
 - b) Measure a Normal Sweep curve in the "PROBE" screen.
 - c) Press F4 and then F7 to enter the "COUPLER" screen.
 - d) Make sure the "MODS" and "RESER" are all 0's.
 - e) Use F6 to set the source from 40-100dB(typically 60dB).
 - f) The Coupler Gain Target can now be viewed. Use DATA/GRAPH to select graph or tabular form.
 - g) The Compression Ratio for each freq can also be viewed if the DSL WDRC rule is used.
- 13) View the Coupler SSPL target:
 - a) Press F5 to enter "SSPL" mode.
 - b) Make sure the "MODS" and "RESERVE" are all 0's.
 - c) The Coupler SSPL Target can now be viewed as both a graph and tabulation.

How to Verify the FP40's Compatability with DSL 4.1

10/31/01

DSL 4.1 Instructions:

- 1) Enter the patients birth date.
 - a) Go to the "Client Information" Window.
 - b) Enter the patients birth date.
- 2) Fill out the Assessment Data Window.
 - a) Go to the "Assessment Data" Window.
 - b) Enter the thresholds into the "Threshold" row.
 - c) If custom UCL is desired, enter them into the "Upper Limit" row.
 - d) If custom RECD is desired, enter them into the "RECD" row.
- 3) Set up the "Assess" Menu:
 - a) "Transducer": Always set to ER3 Insert Earphone.
 - b) "HL to SPL transform":
Set to "Predicted" for average RECD, or to "RECD" for measured RECD.
- 4) Set up the "HA Fit" Menu:
 - a) "Circuit Type": Set to Linear , or WDRC(Fixed CR).
 - b) "Style": Select the aid type(BTE, CIC, ITC, ITE).
Note: CANAL on the FP40 is ITC on DSL 4.1.
 - c) "RE to 2CC transform":
Set to "Predicted"for average RECD, or to "Measured" for measured RECD.
 - d) "Maximum Output":
Set to "Predicted" if using predicted UCL's, or to "Measured" if using measured UCL's.
 - e) "Speech Spectrum":
Set to UWO Child if using any childs age on the FP40.(0-2 M, ..., >5 YR).
Set to Cox and Moore if using "ADULT" age on the FP40.
 - f) "Compression Threshold": Set to (40, 45, 50, 55, 60, or 65).
- 5) View the Real Ear Average Speech-Weighted Input Target:
 - a) Set the "Verify" Menu Mid-Level Input:
 - i) Check both Speech Weighted and Real Ear Aided Response.
 - b) Go to the Verification Mid-Level Targets Window.
 - c) The target can be viewed in the middle row labeled as the selected Speech Spectrum.
- 6) View the Real Ear Maximum Output Target:
 - a) Set the "Verify" Menu High-Level Input to Real Ear Aided Response.
 - b) Go to the Verification High/Low-Level Targets Window.
 - c) The target can be viewed in the top row labeled as "Max. Output".
- 7) View all three Real Ear Targets in graphical form:
 - a) Go to the Speech SPLogram Graph Window.
 - b) Click the Thresholds/Measured, Upper Limits/Measured, and Aided Target/Average boxes.
 - c) The targets can be viewed on the graph.

How to Verify the FP40's Compatability with DSL 4.1

10/31/01

DSL 4.1 Instructions(cont):

- 8) View the Coupler Gain Average Speech-Weighted Input target:
 - a) Set the "Verify" Menu Mid-Level Input:
 - i) Check both Speech Weighted and Coupler Gain.
 - b) Go to the Verification Mid-Level Targets Window.
 - c) The target can be viewed in the middle row labeled as the selected Speech Spectrum.
- 9) View the Coupler Gain Average Constant-Level Input target:
 - a) Set the "Verify" Menu Mid-Level Input:
 - i) Check both Constant-Level and Coupler Gain.
 - b) Go to the Verification Mid-Level Targets Window.
 - c) Set the "Measured" box to the left of the middle row to the same as the FP40 in step 12d.
 - d) The target can be viewed in the middle row.
 - e) This target can also be viewed on the "Hearing Aid Specification" window.
 - i) Set the Res. Gain to 0.
 - ii) Set the Input to the same as the FP40 in step 12d.
 - iii) Read the target on the second and third rows.
 - iv) The compression Ratio for each frequency can be read on the fourth row if WDRC is selected.
- 10) View the Coupler SSPL Target:
 - a) Set the "Verify" Menu High-Level Input to Coupler Output.
 - b) Go to the Verification High/Low-Level Targets Window.
 - c) The target can be viewed in the top row labeled as "Max. Output".
 - d) This target can also be viewed on the "Hearing Aid Specification" window.
 - i) Read the target on the top row labeled as "SSPL".

Ref: FP40_DSL_Verify.doc