


## 8000 BIOS Setup Instructions

1. Turn off the 8000
2. Plug a USB computer keyboard into one of the USB ports on the back of the 8000
3. Press and hold the Delete key on the keyboard
4. Turn on the 8000
5. Continue holding the delete key until the BIOS setup screen appears as shown below



6. Select "Standard CMOS Features" and press enter.
7. Set each item to match the defaults listed on the following pages.
8. Use the up and down arrows to select an item to change then press enter to change it.
9. When the settings match the defaults press ESC.
10. Repeat this process for Advanced BIOS Features, Advanced Chipset Features, Integrated Peripherals, Power Management Setup and PnP/PCI Configurations.
11. When all of the settings match Press F10
12. Press enter to "SAVE to CMOS and EXIT"
13. At this point the 8000 should boot up normally.

# Bios Settings for Fonix 8000

## Standard CMOS Features

Date  
(mm/dd/yyyy):  
Mon Aug 3 2009  
Time (hh/mm/ss):  
08:29:30

IDE Primary  
Master

*[SanDisk  
SDCFH2-  
002G]*

IDE HDD Auto-Detection

[Press  
Enter]

IDE Primary Master  
Access Mode

[Auto]  
[Auto]

IDE Primary  
Slave:

*[None]*

IDE HDD Auto-Detection

[Press  
Enter]

IDE Primary Master  
Access Mode

[Auto]  
[Auto]

Drive A: *[None]*

Drive B: *[None]*

Video *[EGA / VGA]*

Halt on *[All But  
Keyboard]*

## Advanced BIOS Features

Virus Warning	[Disabled]
CPU Internal Cache	[Enabled]
First Boot Device	[USB-ZIP]
Second Boot Device	[HDD-0]
Third Boot Device	[Disabled]
Boot Other Device	[Disabled]
Swap Floppy Drive	[Disabled]
Boot Up Floppy Seek	[Disabled]
Boot Up NumLock Status	[Off]
Gate A20 Option	[Fast]
Typematic Rate Setting	[Disabled]
Typematic Rate (Char/Sec)	6
Typematic Delay (Msec)	250
Security Option	[Setup]
OS Select For DRAM > 64MB	[Non-OS2]
Full Screen LOGO Show	[Disabled]
Small Logo(EPA) Show	[Enabled]

## Advanced Chipset Features

Memory Frequency	[200 MHz]
Video Memory Size	[8M]
Output display	[CRT]
TFT / LVDS Resolution	[800 x 600]
Onboard Audio	[Enabled]
Onboard USB Controller	[Enabled]
Onboard USB2.0	[Enabled]
Onboard USB UDC	[Disabled]
onboard USB OTG	[Disabled]
Onboard IDE	[Enabled]
Memory Hole At 15M-16M	[Disabled]

## Integrated Peripherals

Master Drive PIO Mode	[Auto]
Slave Drive PIO Mode	[Auto]
IDE Primary Master UDMA	[Auto]
IDE Primary Slave UDMA	[Auto]
IDE DMA transfer access	[Enabled]
Onboard Lan1 Control	[Enabled]
Onboard Lan2 Control	[Disabled]
Onboard LAN Boot Rom	[Disabled]
Digital I/O Device	[Press Enter]
IDE HDD Block Mode	[Enabled]
PS/2 Mouse Support	[Enabled]
Onboard FDC Controller	[Disabled]
Onboard Serial Port 1	[Disabled]
Onboard Serial Port 2	[3F8/IRQ4]

Serial Port 2 I/F	[RS232]
UART Mode	[Normal]
Select	
UR2 Duplex	[Half]
Mode	
Onboard Parallel	[378/IRQ7]
Port	
Parallel Port	[SPP]
Mode	
Backlight	[The
controller	Medium
	Light]

## Power Management Setup

ACPI Suspend	[S1(POS)]
Type	
Power	[APM]
Management	
PWRON After	[Off]
PWR-Fail	
MODEMUse IRQ	[N/A]
PME Event	[Enabled]
function	
Soft-Off by PWR-	[Instant-Off]
BTTN	
Power-On by	[Disabled]
Alarm	

## PnP/PCI Configurations

PNP OS Installed	[Yes]
Init Display First	[Onboard]
Reset	
Configuraiton	[Disabled]
Data	
Resources	[Manual]
Controlled By	

## IRQ Resources

[Press  
Enter]

IRQ-3 assigned to	[Reserved]
IRQ-4 assigned to	[PCI Device]
IRQ-5 assigned to	[Reserved]
IRQ-7 assigned to	{PCI Device]
IRQ-10 assigned to	[PCI Device]
IRQ-11 assigned to	[PCI Device]

## Memory Resources PCI/VGA Palette Snoop

[Press  
Enter]

[Disabled]

## ISA Setup

[Press  
Enter]

ISA I/O 3E0- 3EFh	[Enabled]
ISA I/O 2E0- 2EFh	[Enabled]
ISA I/O 100- 13Fh	[Disabled]
ISA I/O 200- 27Fh	[Disabled]
ISA I/O 300- 37Fh	[Disabled]
ISA I/O A79h	[Disabled]
ISA memory CC000- CFFFFh	[Disabled]
ISA memory D0000- D7FFFh	[Disabled]
ISA memory D8000- DFFFFh	[Disabled]