


New FONIX[®] Integrated Probe Microphone

for the 7000 Hearing Aid Test System

Smaller... Lighter... More Adjustable...

- Small detachable reference microphone
- Above the ear placement increases test accuracy
- Raised shelf. More room for the hearing aid!
- Adjustable vertically by 1.5 inches
- Pivots ± 15 degrees for easier probe placement
- Guide for easy probe attachment
- Half the weight of the current probe system!

The new Integrated Probe Microphone comes standard with any new 7000 Hearing Aid Test System with the Real-Ear Option. It can also be retrofitted onto existing 7000 systems.



M707 Integrated probe microphone



FRYE ELECTRONICS, INC.

P.O. Box 23391 • Tigard, OR 97281-3391 • USA

(503) 620-2722 • (800) 547-8209

Fax: (503) 639-0128

www.frye.com • e-mail: sales@frye.com