


Introducing the

FP35 Analyzer External Monitor Option

Frye Electronics, Inc. is pleased to announce the External Monitor Option on the FP35 analyzer. The External Monitor Option allows the user to connect the FP35 analyzer to an external video monitor for easier viewing of the FP35 screen by the user, the client, and any others in the room. When the external monitor is connected, both the monitor and the built-in LCD display are active, letting the user view test results on either screen.

Just about any standard color computer monitor will work, including a modern flat screen. Although, like the built-in LCD display, results are shown in only two colors, the user may choose these colors from a variety of options.

Most existing FP35 analyzers can be upgraded to include this option. Additional hardware upgrades may be required of units with serial numbers below #351.



FRYE ELECTRONICS, INC.

P.O. Box 23391 • Tigard, OR 97281-3391 • USA

(503) 620-2722 • (800) 547-8209

Fax: (503) 639-0128

www.frye.com • e-mail: sales@frye.com