

FONIX[®] FP40 software changes and improvements

<p>Version 1.0 (Original Release 3/90)</p>	<p>Version 3.1 (1/95) Added: AI (Audibility Index). New Laser Printout format. External Color Ink Jet Printer Interface.</p>
<p>Version 1.2 (5/90) Added: JIS Standard.</p>	<p>Version 3.2 (4/95) Added: CIC and ANSI '92</p>
<p>Version 1.3 (8/90) Added: ISI Standard; CRT Option.</p>	<p>Version 3.23 (4/96) Added: #5 battery pill as optional accessory; label for #5 battery pill.</p>
<p>Version 2.0 (3/91) Added: Fast Sweep</p>	<p>Version 3.30 (11/96) Added: Real ear total SPL output. New sound chamber and probe speaker.</p>
<p>Version 2.1 (2/92) Added: New Multi-Curve in Real Ear & Coupler. Target 2cc FOG & SSPL90. Easy-to-read Double-Size Characters. Full & Partial Menu Choices. Probe Signal Selection (F7 Definition).</p>	<p>Version 3.40 (6/97) Added: DSL 4.1 Linear and WDRC.</p>
<p>Version 2.3 (4/93) Added: 40 dB Input Source (K-Amp). Spectrum Analysis Mode. Improved: Video Display Format (everything on one screen). Target 2cc FOG vs Actual FOG Comparison in Test Box. Equivalent Input Noise (EIN) Measurement.</p>	<p>Version 3.50 (10/99) Added: Digital speech testing and expanded DSL</p>
<p>Version 2.4 (10/93) Added: Screen saver</p>	<p>Version 3.60 (4/00) Added: ANSI '96 and Profiler automated test sequences.</p>
<p>Version 2.5 (2/94) (2.6 Last ver. available for old CPU boards) Added: RS232 Computer Interface. Laser Printer Interface. CHAP & HearWare Interface.</p>	<p>Version 3.70 (9/05) Added: ANSI '03</p>
<p>Version 3.02 (10/94) Added: New CPU, VGA Color Option, Telecoil Option, (Earlier FP40/FP40-D can be upgraded.) Interfaces with CHAP and HearCare, a new office management software package. Stores right/left audiograms and targets.</p>	<div data-bbox="846 1761 990 1904" data-label="Image"> </div> <p>FRYE ELECTRONICS, INC. P.O. Box 23391 • Tigard, OR 97281-3391 • USA (503) 620-2722 • (800) 547-8209 Fax: (503) 639-0128 www.frye.com • e-mail: sales@frye.com</p>