

FONIX Price List

U.S. and International Prices



FRYE ELECTRONICS, INC.
 P.O. Box 23391 • Tigard, OR 97281-3391 • USA
 (503) 620-2722 • (800) 547-8209
 Fax: (503) 639-0128
 www.frye.com • e-mail: sales@frye.com

October 1, 2011

See our online Price List for up-to-date prices and pictures of FONIX instruments and accessories.

F.O.B. TIGARD

Prices subject to change without notice • \$25.00 Minimum order • 15% Restocking fee

Windows software programs (such as the NOAH Real-ear Module and WinCHAP), and audiometric CDs (such as ANL and the Audiometric Test CD) are not returnable.
 Warranty is one year parts and labor on new instruments, 6 months on used instruments and repairs.

FONIX 8000 HEARING AID TEST SYSTEM

FONIX 8000 Module	8000	\$10,000
Includes one of the following test sequences: ANSI 03/96/87 (including Telewand), IEC 05/94, or JIS.		
<i>Note: Monitor sold separately</i>		

Chamber Packages with FONIX 8000 module purchase

8120 polar plot chamber	FC8120	\$4,949
Includes M1958E coupler microphone, magnetic HA-1, HA-2, and CIC couplers, ear level adapter, auto detecting battery simulator pills (#13, #675, #10 and #312), rotator accessories for BTE and ITE directional tests, blue stik, extra paper, USB computer cable, and removable test surface.		
8050 chamber (no Polar Plot)	FC8050	\$3,600
Includes M1958E coupler microphone, HA-1, HA-2, and CIC couplers, ear level adapter, auto detecting battery simulator pills (#13, #675, #10 and #312), blue stik, extra paper, and USB computer cable.		

Real-Ear for FONIX 8000 HEARING AID TEST SYSTEM

Probe microphone (M767) with an advanced integrated ear hook design, sound field speaker, either a floor stand or a swing arm (which can be attached to a wall or desk), Visible Speech, Improved Digital Speech signal for testing open fit hearing aids with probe remote module, monitor headset, calibration adapters, probe tubes, and marking pen.	8000 PROBE/STAND OR 8000 PROBE/ARM	\$2,700
Same as above with additional M767 microphone for Dual probe. Special price when purchased with 8000	8000 PROBE/STAND/DUAL OR 8000 PROBE/ARM/DUAL	\$3,450

Options ordered with purchase

Additional test sequences (Specify ANSI (including telewand), IEC, or JIS)		\$400 each
OES Option including MZ couplers	8000 OES	\$500
OES Option (No MZ couplers)	8000 OES w/o	\$100
VA-CORFIG	8000 VA-CORFIG	\$400

Accessories

Computer Programs (<i>See the Computer Programs price sheet for more details</i>)		
FONIX NOAH Real-ear Module	FNMV2.21AN	\$400
WINCHAP	WINCHAP V3.00	\$600
FONIX Troubleshooter	TSHOOT V1.11	\$400
Android Remote Keyboard: Includes Android tablet and pre-installed FONIX 8000 Keyboard App. Requires 8000 software v2.14.	8000 BLUEKEY	\$495
RECD package. Includes insert earphone with assortment of eartips.	015-3418-01	\$425
LCD Monitor 17"	058-0078-00	\$550
External Telecoil board and cable	120-1005-00	\$275

Hearing instrument programming cable (HiPro cable), 6ft.	072-0317-00	\$27
Speaker stand, tabletop model	030-0021-00	\$50
Demonstration Ear with stand	180-0018-00	\$350
Auto Detect Battery Simulator pill set (includes #10, #13, #675, #312)	180-0022-00	\$240
Chamber Remote Module (for connection to sound chambers other than the 8050/8120)	034-1263-00	\$600
Additional Accessories		
Infant/Child Headband Package including infant, child, and adult headbands, six bendable earhooks, and two sets of "animal ears."	180-0017-01	\$125
Consumables		
Paper roll	026-0006-00	\$8
Probe tubes, <i>package of 25</i>	074-0020-00	\$44
Insert earphone accessories		
Package of 50 medium eartips	025-0670-00	\$35
Package of 50 small eartips	025-0671-00	\$35
Package of 24 large eartips	025-0672-00	\$35
Sound tube replacement	042-0526-00	\$12
Nipple for sound tube	042-0532-00	\$5
Dry erase marker	006-0065-00	\$3
Blue Adhesive Putty	033-0002-00	\$4
Battery Simulator Pills (12" cable), <i>Autodetecting</i>		
#675, <i>Blue</i>	059-3004-00	\$65
#13, <i>Orange</i>	059-3005-00	\$65
#312, <i>Brown</i>	059-3006-00	\$65
#10, <i>Yellow</i>	059-3010-00	\$65
#AA, <i>Black</i>	059-3011-00	\$130
Add Options Post Manufacture		
Additional test sequence, ANSI (including telewand), IEC, or JIS		\$500 each
OES with MZ-1, MZ-2, and MZ-3 couplers	8000 OES/UP	\$750
OES Option (No MZ couplers)	8000 OES w/o/UP	\$250
VA-CORFIG	8000 VA-CORFIG/UP	\$500

Replacement Parts		
<i>Real-ear Equipment</i>		
Probe Remote Module	034-2263-00	\$800
Speaker swing arm package	ACC/SWING/8000	\$495
Real-ear speaker	034-2261-00	\$250
Speaker swing arm without speaker. Includes cable with 3.5 mm plug	034-4456-01	\$345
Monitor headset with 1/4" adapter	015-1214-00	\$47
Speaker stand, floor model	030-0020-02	\$150
Phone plug adapter, 1/8" to 1/4", for monitor headset	044-0042-00	\$6
<i>Sound chamber Equipment</i>		
8120 Sound chamber	FC8120	\$4,949
8050 Sound chamber	FC8050	\$3,600
Telewand for ANSI 96/03 telecoil measurements	043-1053-00	\$115
Sound Chamber Cable (8050/8120)	119-0282-00	\$100
Clamp Assembly BTE (8120)	034-1974-00	\$75
BTE Clamp Liners (set of 5) (8120)	011-0002-00	\$6
Rotator Post ITE (8120)	053-1012-00	\$75
Rotator Post BTE (8120)	053-1013-00	\$75
Removable Test Surface (8120)	053-1014-00	\$70
<i>Couplers & Adapters</i>		
HA-1 2-cc direct access coupler	044-1031-01	\$75
HA-2 2-cc coupler	044-1029-01	\$83
HA-1 Magnetic Coupler (for the 8120)	044-5031-01	\$75
HA-2 Magnetic Coupler (for the 8120)	044-5029-01	\$83
Adapter 14 mm to 1 in (for sound level calibrators)	044-1006-01	\$30
Ear Level Adapter	044-1003-00	\$20
CIC coupler for CIC Option	044-1027-01	\$75
Open Fit coupler	044-1023-01	\$100
Rectangular reference microphone adapter	044-0208-00	\$45
Probe adapter, 14mm to probe	044-0025-01	\$25
<i>OES Option Couplers</i>		
MZ-1 coupler	044-1024-01	\$82
MZ-2 coupler	044-1022-01	\$95
MZ-3 coupler	044-1021-01	\$82
Ear level adapter for HA and MZ couplers	044-1003-00	\$20
<i>Misc Couplers and Related Repair Parts</i>		
HA-3 modified coupler	044-1015-04	\$120
SB-1 2-cc for universal-fit hearing aids, 9-11mm dia.	044-2031-01	\$95
6-cc coupler, NBS 9A for testing earphones	044-1004-03	\$250
HA-1 coupler for B&K 1" microphone	044-1009-00	\$500
Adapter sleeve for B&K HA-1 coupler	044-1008-00	\$100
Ear Hook assembly for M767 Microphone		
Slide Ear Hook	044-0070-00	\$35
Ear Hook with Velcro	043-2051-04	\$5
#13 tubing for ear level adapter, price per inch	074-0006-00	\$1

"O" Ring set of 3 for HA-2	042-1500-00	\$3
"O" Ring set of 5 large	042-1502-00	\$5
"O" Ring set of 4 for all couplers	042-1504-00	\$6
Microphones		
M1958E Elec. Microphone with 18in. cable	M1958E	\$750
M1950E Elec. Microphone with 6ft. cable	M1950E	\$750
M767 Integrated probe microphone	M767	\$1,100
Dummy Microphone	041-0025-00	\$9
Other		
8000 Operator's Manual	010-0035-01	\$45
8000 Maintenance Manual	010-0034-00	\$45
Cable 6' USB ver.2 type a/b for computer connection	072-0313-00	\$15

FONIX Price List

U.S. and International Prices



FRYE ELECTRONICS, INC.
 P.O. Box 23391 • Tigard, OR 97281-3391 • USA
 (503) 620-2722 • (800) 547-8209
 Fax: (503) 639-0128
 www.frye.com • e-mail: sales@frye.com

October 1, 2011^A

See our online Price List for up-to-date prices and pictures of FONIX instruments and accessories.

F.O.B. TIGARD

Prices subject to change without notice • \$25.00 Minimum order • 15% Restocking fee

Windows software programs (such as the NOAH Real-ear Module and WinCHAP), and audiometric CDs (such as ANL and the Audiometric Test CD) are not returnable.
 Warranty is one year parts and labor on new instruments, 6 months on used instruments and repairs.

FP35 Portable Hearing Aid Analyzer



FONIX FP35 Hearing Aid Analyzer

FONIX FP35 Coupler Package

Software Includes: pure-tone coupler multicurve, composite and one of the following test sequences: ANSI 03, ANSI 96, IEC 05, IEC 94, JIS, or ISI.

Hardware includes: coupler microphone, HA-1 and HA-2 couplers, calibration adapter, ear-level adapter, Blue Stik, RS232 with serial to USB cable, thermal printer with roll of paper, line-in connector, built-in display, external parallel printer port, with medical grade power supply.
 Battery Current Drain: Battery set sold separately.

FP35 MG	\$4,785
---------	---------

FONIX FP35 with Real-ear

Coupler package plus **Software Includes:** the Real-ear Option feature for probe microphone measurements including Visible Speech.

Hardware includes: integrated probe microphone(M767), probe tubes, marking pen, calibration clip, rectangular reference microphone adapter, and speaker post with extension post, with medical grade power supply.
 Battery Current Drain: Battery set sold separately.

FP35/P MG	\$5,620
-----------	---------

Options ordered with purchase

Additional test sequence (Specify ANSI 96, ANSI 03, JIS, IEC 05, IEC 94, or ISI)		\$125 each
CIC Option including coupler	FP35 CIC	\$100
OES including MZ couplers	FP35 OES	\$500
OES Option (No MZ couplers)	FP35 OES w/o	\$100
External Monitor Option (not including monitor)	FP35 EXT. Monitor	\$500
Screening Audiometer Option including software and earphone calibration adapter	FP35 AUD	\$150
Screening Audiometer Option plus ER3A earphones	FP35 AUD ER3A/08	\$625
Spanish/French Language	FP35 V6.30S/F	no charge

Accessories

Computer Programs (See the Computer Programs price sheet for more details)

FONIX NOAH Real-ear Module	FNMV2.21AN	\$400
WINCHAP	WINCHAP V3.00	\$600
FONIX Press & Go	PRESS&GO	\$300
FONIX Troubleshooter	TSHOOT V1.11	\$400
LCD Monitor 17"	058-0078-00	\$550
Auto Detect Battery Simulator pill set (includes #10, #13, #675, #312)	180-0022-00	\$240
Auto Detect Battery Simulator pill #AA, Black	059-3011-00	\$130
iPod with pre-loaded Audiometric Speech Tests. Includes docking station and cable for connection to the FP35	015-2010-01	\$299
Soft Carrying Case	063-0007-00	\$189

Telewand for ANSI 96/03 telecoil measurements	043-1050-00	\$115
Demonstration Ear with stand	180-0018-00	\$350
Additional Accessories		
Ear Hooks		
Wedge style	043-2043-00	\$40
Children's wedge style	043-3043-00	\$40
Infant/Child Headband Package including infant, child, and adult headbands, six bendable earhooks, and two sets of "animal ears."	180-0017-01	\$125
Old style (no wedge)	034-1776-00	\$28
Children's old style (no wedge)	034-1777-00	\$28
Headband	042-0516-00	\$40
Additional Accessories for FP35 analyzers (board version 08 and above), Serial No. required		
Speaker Swing Arm Package for FP35 with 3.5 mm plug. Includes speaker, swing arm, cable, and mounting brackets.	ACC/SWING/ FP35/08	\$495
Speaker Stand Package for FP35 with 3.5 mm plug. Includes speaker, speaker stand, and cable.	ACC/FP35/SPE/08	\$400
Internal Speaker Cable (1ft long, stereo 3.5mm both ends)	119-0363-00	\$9
RECD package for 1/4 in plug. Includes ER3 insert earphone with assortment of eartips	015-4418-00	\$450
External sound chamber, 6040. Requires M311 microphone with a six foot cable if ordered with a new FP35 analyzer	FC6040/FP35/08	\$1,500
External sound chamber, 6050 with 3.5 mm battery pill jack. Requires M311 microphone with a six foot cable if ordered with a new FP35 analyzer	FC6050/FP35/08	\$1,500
Telecoil board with 3.5 mm plug	120-1003-00	\$200
Cable, Y 1/4" stereo plug to 2 1/4" mono plugs (RECD insert Y cable package)	119-0344-00	\$10
Screening Audiometer Option including ER3A earphones	FP35 AUD ER3A/08	\$625
Screening Audiometer Option with software and earphone calibration adapter	FP35 AUD	\$150
Cable, FP35 to 6050/7020 sound chamber	119-0186-00	\$50
Additional Accessories for FP35 analyzers (board version 07 and below), Serial No. required		
Speaker Swing Arm Package for FP35 with modular plug. Includes speaker, swing arm, cable, and mounting brackets.	ACC/SWING/FP35	\$495
Speaker Stand Package for FP35 with modular plug. Includes speaker, speaker stand, and cable.	ACC/FP35/SPE	\$400
Internal Speaker Cable (1ft long, 4pin modular both ends)	072-0254-00	\$4
RECD package for 3.5 mm plug. Includes ER3 insert earphone with assortment of eartips.	015-2418-01	\$440
External sound chamber, 6040 with modular plug	FC6040/FP35	\$1,500
Y-Cable, 3.5 mm x 1/4 " module plug for Dual Insert Earphones	119-0341-00	\$60
Screening Audiometer Option including ER3A earphones	FP35 AUD ER3A	\$625
Screening Audiometer Option with software and earphone calibration adapter	FP35 AUD	\$150
Consumables		
Paper roll	026-0019-00	\$7
Probe tubes, <i>package of 25</i>	074-0020-00	\$44
Insert earphone accessories		
Package of 50 medium eartips	025-0670-00	\$35
Package of 50 small eartips	025-0671-00	\$35
Package of 24 large eartips	025-0672-00	\$35
Sound tube replacement	042-0526-00	\$12
Nipple for sound tube	042-0532-00	\$5

Dry erase marker	006-0065-00	\$3
Blue Stik (previously called "Fun Tak")	033-0002-00	\$4
Battery Simulator Pills (12" cable), <i>Autodetecting</i>		
#675, <i>Blue</i>	059-3004-00	\$65
#13, <i>Orange</i>	059-3005-00	\$65
#312, <i>Brown</i>	059-3006-00	\$65
#10, <i>Yellow</i>	059-3010-00	\$65
#AA, <i>Black</i>	059-3011-00	\$130
Software Upgrade/Add Options Post Manufacture, <i>Serial No. required</i>		
OES, including MZ couplers	FP35 OES	\$500
Composite/Digital Speech	FP35 COMP/UP	\$750
Additional test sequence, ANSI 96, ANSI 03, IEC 94, IEC 05, or JIS		\$175 each
Spanish/French Language (upgrade from English to Spanish/French)	FP35 V6.30S/F	\$300
Software upgrade (within the same languages)	FP35 V6.30E(S)	\$200
Software upgrade with new hardware	Contact Frye for information.	
Hardware Upgrades (<i>This upgrade must be performed at Frye factory.</i>)		
Include Board-10, new microphone board, and modify existing enclosure. <i>Mic. sold separately</i>	FP35UP	\$1,750
Upgrade from coupler to real ear, includes software, mic (new style) (<i>for units with serial numbers above 2000</i>)	FP35PROBE/UP2	\$1,100
Upgrade from coupler to real-ear, includes new M300 & requires FP35UP (<i>for units with serial numbers under 2000</i>)	FP35PROBE/UP	\$1,100
Replacement parts		
<i>Real-ear equipment</i>		
External Real-ear speaker	034-2261-00	\$250
Probe speaker cables FP35 analyzers (board version 08 and above), 3.5 mm plug		
10 ft	119-0176-02	\$50
25 ft	119-0187-01	\$50
50 ft	119-0185-00	\$75
Swing arm and cable with 3.5 mm plug without speaker (board version 08 and above)	034-4456-01	\$345
Probe speaker cable for FP35 (board version 07 and below), modular plug		
10 ft	119-0179-01	\$50
50 ft	119-0180-01	\$60
Swing arm and cable with modular plug without speaker (board version 07 and below)	034-2456-01	\$345
Speaker post extension	034-1533-02	\$25
Mounting sleeve for probe microphone	044-1204-00	\$15
Mounting sleeve for reference microphone	044-1206-00	\$15
Probe calibration clip	044-0200-00	\$18
Speaker stand, floor model	030-0020-02	\$150
1/4" to 1/8" adapter for use with insert earphones	044-0042-00	\$5
<i>Standard Couplers</i>		
HA-1 2-cc direct access coupler	044-1031-01	\$75
HA-2 2-cc coupler	044-1029-01	\$83
CIC coupler for CIC Option	044-1027-01	\$75
Open Fit coupler	044-1023-01	\$100
Rectangular reference microphone adapter	044-0208-00	\$45

Calibrator adapter, 14mm to 1 in	044-1006-01	\$30
Probe adapter, 14mm to probe	044-0025-01	\$25
Probe calibration clip	044-0200-00	\$18
OES Option Couplers		
MZ-1 coupler	044-1024-01	\$82
MZ-2 coupler	044-1022-01	\$95
MZ-3 coupler	044-1021-01	\$82
Ear level adapter for HA and MZ couplers	044-1003-00	\$20
Misc Couplers and Related Repair Parts		
Earphone calibration adapter	044-1020-01	\$30
SB-1 2-cc for universal-fit hearing aids (9-11mm dia.)	044-2031-01	\$95
#13 tubing for ear level adapter (price per inch)	074-0006-00	\$1
“O” Ring set of 3 for HA-2	042-1500-00	\$3
“O” Ring set of 5 large	042-1502-00	\$5
“O” Ring set of 4 for all couplers	042-1504-00	\$6
Microphones		
M767 Integrated probe & reference microphone set	M767	\$1,100
M311 coupler microphone with a 6 ft cable. For use with external sound chambers	M311	\$400
M367 coupler microphone	M367	\$400
Other		
Replacement strap for soft carrying case	063-0008-00	\$15
USB serial cable, 5 meters	072-1311-00	\$90
RS232 cable (integrated) for '08 board design only	119-1400-00	\$20
RS232 cable kit with LED connector (07 Board)	119-0003-00	\$75
Printer cable	119-0326-01	\$27
FP35 Operator's Manual	010-0013-10	\$30
FP35 Spanish Operator's Manual	010-0027-04	\$30
FP35 Maintenance Manual	010-0017-00	\$30
Calibration		
ISO calibration	FP35 CAL/CE	\$225

FONIX Price List

U.S. and International Prices



FRYE ELECTRONICS, INC.
 P.O. Box 23391 • Tigard, OR 97281-3391 • USA
 (503) 620-2722 • (800) 547-8209
 Fax: (503) 639-0128
 www.frye.com • e-mail: sales@frye.com

October 1, 2011

See our online Price List for up-to-date prices and pictures of FONIX instruments and accessories.

F.O.B. TIGARD

Prices subject to change without notice • \$25.00 Minimum order • 15% Restocking fee

Windows software programs (such as the NOAH Real-ear Module and WinCHAP), and audiometric CDs (such as ANL and the Audiometric Test CD) are not returnable. Warranty is one year parts and labor on new instruments, 6 months on used instruments and repairs.

FONIX Audiometers

FA10 Digital Hearing Evaluator

Including 110 dB output, Output Reverse Master Hearing Aid Simulator, Stenger test, B71 bone vibrator headset, Audiogram pad, SISI, ABLB, and MLB



FONIX FA-10 shown with optional gooseneck microphones

With TDH-39 Earphone Headset	FA10 STD	\$3,495
With 3A Insert Earphones	FA10/INS	\$3,595
Dual Calibration with TDH-39 and 3A inserts	FA10/DUAL	\$3,995

Audiometer RS232 upgrade option (in house). Includes internal instrument hardware and external USB cable for RS232 communication. Communication to a Frye audiometer requires the RS232 option be installed on the instrument. This is independent of the software purchase such as Winchap or Fonix Noah Module.

FA RS232 add \$250

FA12 Digital Hearing Evaluator

Including 120 dB output, Master Hearing Aid Simulator, Stenger test, B71 bone vibrator headset, Audiogram pad, SISI, ABLB, and MLB

With TDH-39 Earphone Headset	FA12/STD	\$3,495
With 3A Insert Earphones	FA12/INS	\$3,595
Dual Calibration with TDH-39 and 3A inserts	FA12/DUAL	\$3,995
Audiometer RS232 upgrade option (in house). Includes internal instrument hardware and external USB cable for RS232 communication. Communication to a Frye audiometer requires the RS232 option be installed on the instrument. This is independent of the software purchase such as Winchap or Fonix Noah Module.	FA RS232	add \$250

Accessories Packages

Clinical Office Package

Including Boom microphone & monitor headset, Talkback microphone, Desktop stand and clip for talkback microphone, Dust Cover, Patient Response Switch, 2 sound field speaker & 10 ft cables, 2 speaker wall mounts, 7 6' length patch cords 1/4" to 1/4" plug.

Special Package Price (save \$244)

FA CLINIC PKG \$925

Sound booth package

Including Audio headset adapter, Patient Response Switch, Economy Talkback microphone with adjustable stand, Monitor headset with adapter, 5 6' length patch cords 1/4" to 1/4" plug.

Special Package Price (save \$32)

180-0100-00 \$250

Accessories

Computer Programs (See the Computer Programs price sheet for more details)

FONIX NOAH Audiometer Module	FNMV2.02AU	\$400
WINCHAP	WINCHAP V3.00	\$600

FONIX AUDIOMETERS [10]

Boom microphone & monitor headset	015-2217-01	\$190
Talkback microphone	058-1027-00	\$75
Desktop stand and clip for talkback microphone	030-1026-00	\$60
Dust Cover	063-0056-00	\$18
Patient Response Switch	076-1700-02	\$60
Soundfield speaker package (2 sound field speakers & 10 ft cables)	041-2503-00	\$475
2 speaker wall mounts	030-1028-00	\$95
Soft carrying case	063-1011-00	\$195
Acceptable Noise Level (ANL) Test, CD/Audiometer	ANLTESTCD	\$55
Audiometric Speech Tests CD and word list	027-0023-02	\$50
TEN (HL) Audiometric Tests CD	TEN HL	\$40
iPod with pre-loaded Audiometric Speech Tests. Includes docking station and y-cable for connection to the audiometer	015-1010-01	\$299
CD Velcro mounting kit	034-2912-01	\$20
Pediatric headphone headband, Wilfan,	015-0222-00	\$80
Monitor headset with adapter	015-1214-00	\$47
Gooseneck microphone package including 2 gooseneck microphones and monitor headset	GOOSENECKOPT	\$235
Economy Talkback microphone with adjustable stand	058-1015-00	\$25
Telescopic stand and clip for talkback microphone	030-1021-00	\$90
Word list, 50 sheets, 2 sheets per set	010-0054-01	\$11
Upgrades		
Audiometer RS232 field upgrade option. Includes internal instrument hardware and external USB cable for RS232 communication. Communication to a Frye audiometer requires the RS232 option be installed on the instrument. This is independent of the software purchase such as Winchap or Fonix Noah Module.	FA RS232 KIT	\$250
Dual calibration kit including software, operator's manual, and label	FA DUAL/KIT	\$125
In-house upgrade to dual calibration including necessary earphones calibration, software, and label		
Insert Earphones	FA UPG DUAL ER3A	\$875
TDH Headphones	FA UPG DUAL TDH	\$775

Replacement Parts

Main Components

B71 Bone oscillator

Complete assembly	015-1224-00	\$200
B71 100 ohm bone vibrator	015-0226-00	\$135
B71 cord	015-0412-00	\$50
B71 headband	015-0224-00	\$30
3A 50 ohm insert earphones		
Complete assembly including eartips	015-0416-00	\$575
Package of 50 medium eartips	025-0670-00	\$35
Package of 50 small eartips	025-0671-00	\$35
Package of 24 large eartips	025-0672-00	\$35
Sound tube replacement	042-0526-00	\$12
Nipple for sound tube	042-0532-00	\$5
Cord for 3A insert earphones	072-0301-00	\$100
Clip Velcro button	015-0417-00	\$10

TDH-39 Headphones

Complete assembly	015-1221-00	\$525
TDH-39 100 ohm earphone	015-0206-00	\$135
Headphone headband	015-0221-00	\$90
Headphone cord	015-0410-01	\$75
Wilfan Headphones, Complete assembly	015-1222-00	\$425

Cables & Adapters

1/8" stereo to 1/4" stereo plug (<i>Audio headset adapter</i>)	044-0026-00	\$10
Y-cable adapter, 3.5 mm to dual 1/4" plugs for CD player	119-0272-01	\$25
Patch cord cable, 1/4" to 1/4", 6' length	072-0268-00	\$28
Speaker cable, RCA to 1/4" plug, 6' length	072-0260-00	\$18
Speaker cable, RCA to 1/4" plug, 10' length	072-0262-00	\$25
USB serial cable, 5 meters (<i>RS232 options required</i>)	072-1311-00	\$90
RS232 cable kit with LED connector	119-0003-00	\$75

Microphones

Gooseneck microphone, 2 required	058-1026-01	\$100 each
Talkback microphone stand clip (<i>only holder, no stand included</i>)	058-0034-00	\$25

Other

Sound field speaker, 8 ohm, 2 required	034-1261-00	\$250 each
Speaker wall mounts, pair	030-1028-00	\$95
Audio cups, special order	015-0223-00	\$400
Audiogram pad, 50 sheets	010-0052-00	\$12
Carrying case replacement strap	063-0012-00	\$15
FA10/12 Operator's manual	010-0084-00/010-0086-00	\$30
FA10/12 Maintenance manual	010-0001-00/010-0008-00	\$30

Calibration

Dual calibration	FA-XX CAL DUAL	\$400
Single calibration	FA-XX CAL SINGLE	\$300

FONIX Price List

U.S. and International Prices



FRYE ELECTRONICS, INC.

P.O. Box 23391 • Tigard, OR 97281-3391 • USA

(503) 620-2722 • (800) 547-8209

Fax: (503) 639-0128

www.frye.com • e-mail: sales@frye.com

August 1, 2011

See our online Price List for up-to-date prices and pictures of FONIX instruments and accessories.

F.O.B. TIGARD

Prices subject to change without notice • \$25.00 Minimum order

Windows software programs (such as the NOAH Real-ear Module and WinCHAP), and audiometric CDs (such as ANL and the Audiometric Test CD) are not returnable. Warranty is one year parts and labor on new instruments, 6 months on used instruments and repairs.

Computer Programs

All FONIX computer programs require RS232 functionality. This is standard on the 8000/7000/FP35, but is an optional feature that must be purchased on the 6500-CX/FP40 hearing aid analyzers and all FONIX audiometers. The FONIX equipment must also have the following software version:

- 8000: all versions
- 7000: 1.20 or higher
- 6500-CX: 4.60 or higher
- FP40: 3.50 or higher
- FP35: 2.0 or higher
- FA10/FA12: all versions

We recommend the most current software version for best functionality. The computer used for connecting the FONIX equipment must have a free serial port or use a USB to serial adapter. There have been connection issues with computers using Bluetooth technology.

WinCHAP computer hearing aid program

1-2 computers	WINCHAP V3.00	\$600
3-5 computers	WINCHAP SITE	\$900
WinCHAP upgrade to the latest major version	WINCHAP V3.00UP	\$100
FONIX Press & Go, 1-2 computers	PRESS&GO V.1.10	\$300
FONIX NOAH Real-ear Module		
1-2 computers	FNMV2.21AN	\$400
3-5 computers	FNMV2.21 SITEAN	\$550
FONIX NOAH Audiometer Module		
1-2 computers	FNMV2.02AU	\$400
3-5 computers	FNMV2.02SITEAU	\$550
FONIX NOAH Module upgrade to latest major version	FNMV2.XXAXUP	\$100
FONIX Troubleshooter	TSHOOT V1.11	\$400
FONIX Troubleshooter upgrade to latest major version	TSHOOTUP V1.11	\$100

Replacement Cables

RS232 cable (integrated) for 7000/FP35-08 board	119-1400-00	\$20
RS232 cable kit with LED connector for FP40/FA10/FA12/FP35-07 and older	119-0003-00	\$75
USB serial cable, 5 meters	072-1311-00	\$90
USB serial cable for 6500, 5 meters kit	072-2311-00	\$90