


Wedge-Style Earhooks



Standard size
(043-2043-00) \$40

Child's size.
(043-3043-00) \$40

These handy accessories allow for the placement of the reference microphone and the probe microphone on one piece. It makes the use of the headband unnecessary. The small size ear hook was recently added to the optional accessories of the FONIX FP40/FP40-D and FONIX 6500 series with Quik-Probe.

 **FRYE
ELECTRONICS, INC.®**

P.O. Box 23391 • Tigard, Oregon 97281-3391 • USA
Phone (503) 620-2722 • (800) 547-8209 • FAX (503) 639-0128