

8050 Mid-Sized Sound Chamber



The 8050 Mid-Sized Sound Chamber is available for the FONIX 8000 Hearing Aid Test System. It is compact and attractive yet large enough to provide good sound isolation with minimal sound reflections. The speaker is mounted vertically at the back of the chamber, making it easy to position the hearing aid with respect to the speaker. All cables required for normal testing with the 8000 Hearing Aid Test System including the microphone, battery pill, and a Hi-Pro pass-through are connected inside the chamber, so there are no leaks in the seals when the lid is closed.

Note: The Mid-Sized Chamber does not have polar plot measurement capability.

Exterior Size: 11.5"W x 13.5"D x 9.2"H

Weight: 21.8 lbs

Test Area Size: 5"W x 6"D x 1.5"H

Ambient Noise Isolation: 45 dB at 1 kHz



FRYE ELECTRONICS, INC.

P.O. Box 23391 • Tigard, OR 97281-3391 • USA

(503) 620-2722 • (800) 547-8209

Fax: (503) 639-0128

www.frye.com • email: sales@frye.com