

The Unique FP35 *Touch*



Experience Frye Electronics new FP35 Touch! It brings the world of interactive touch panels to hearing aid analysis—a whole new, more intuitive operator interface.

Press the screen to get a complete roadmap; *Touch* the desired destination to navigate there.

Touch Source SPL to set a new level.

Touch “left” (ear) to change the ear.

Swipe down to show a menu.

Touch to choose a menu item.

Swipe up to close the menu.

Plus many more options.... See the on-line operator’s manual to get a complete listing of “Touch” commands.

Interactive FP35 touch panel processing is similar to the popular Android and Apple products. The user thus experiences a whole new instrument interface, making the FP35 analyzer unique in the hearing industry.



FRYE ELECTRONICS, INC.

P.O. Box 23391 • Tigard, OR 97281-3391 • USA

(503) 620-2722 • (800) 547-8209

Fax: (503) 639-0128

www.frye.com • email: sales@frye.com