

Bill of Materials

Assembly		Description	Rev		Stock UOM				
HEARLAB-D		HEARLAB / COMPLETE W/DESKTOP COMPUTER			EA				
Item	Component and Description		Rev	Qty Per	UOM	Obs. Date	Location Code	ISC	Order Quantity
	015-0224-00	HEADBAND/BONE VIBRATOR/ AUDIOMETER		1.00	ea		Aisle8	B	100.00
	015-0226-00	VIBRATOR/ BONE/ AUDIOMETER, B71, SERIALIZED		1.00	ea		Aisle2	B	100.00
	015-0412-00	CORD/ BONE VIBRATOR/ AUDIOMETER Ø		1.00	ea		Aisle2	B	100.00
	<u>Component Part Notes</u> Make sure vendor know to twist tie individually on the next order. RRD 15MAY2009								
	015-0416-00	EARPHONES/ INSERT/ ER-3A (50 OHMS)		1.00	ea		Aisle2	B	25.00
	034-1261-00	SPEAKER ASSEMBLY, SOUNDFIELD LOUDSPEAKER		1.00	ea		Aisle10	M	1.00
	042-0011-00	FEET/ ADHESIVE		4.00	EA		Aisle3	B	1,000.00
	044-0047-00	ADAPTER/ DUAL BANANA PLUG TO RCA PHONO JAC		1.00	ea		Aisle3	B	50.00
	061-0194-00	ECB WIRED HEARLAB IMPEDENCE TEST JIG		1.00	ea		Aisle4	M	1.00
	063-0013-00	ENCLOSURE/PACTEC/SM HEARLAB/ CL-400		1.00	EA		Aisle10	B	10.00
	072-0262-00	CORD/ AUDIO PATCH/ 1/4 IN. PHONE TO RCA/ 10 FT		1.00	ea		Aisle5	B	100.00
	072-0289-00	CABLE, S-VIDEO 12FT		1.00	ea			B	1.00
	072-0312-00	CABLE,10' FIREWIRE 6 PIN HEARLAB IEE1394		1.00	ea		Aisle6	B	1.00
	072-0313-00	CABLE, 6' USB VER.2 TYPE A/B HEARLAB		1.00	ea		Aisle6	B	1.00
	ACC/HEARLAB	ACCESSORIES/ HEARLAB		1.00	EA		SHIP ASSY	M	1.00
	ACTELEC	ELECTRODE / ACTIVE / MOLDED / HEARLAB		1.00	EA		SHIP SWALL	M	100.00
	DESKTOP	COMPUTER/DESKTOP/TESTED OPTIPLEX 330		1.00	EA		test	M	1.00
	ELECPROC	ELECTRODE PROCESSOR / HEARLAB		1.00	EA		assembly	M	10.00
	<u>Component Part Notes</u> This is the Electrode Processor. Two boards (originally called EPU upper and lower). All electrodes plug into this module. This module attaches to the Stimulus Controller.								
	GNDELEC	ELECTRODE / GROUND / MOLDED / HEARLAB		1.00	EA		SHIP SWALL	M	100.00
	REFELEC	ELECTRODE / REFERENCE / MOLDED / HEARLAB		1.00	EA		SHIP SWALL	M	100.00
	REFMIC	MICROPHONE/ REFERENCE / HEARLAB		1.00	EA		test	M	1.00
	STIMCONT	STIMULUS CONTROLLER (MAIN HEARLAB CHASSIS)		1.00	EA		assembly	M	10.00
	<u>Component Part Notes</u> This is the main chassis. Two large boards (originally called SRU upper and lower) plus the new mepos power supply.								