

Document Maintainer: Kristina Frye

Release Date Document Change Log

November 19, 2009 Initial release

1.0 Purpose

The purpose of this document is to describe the physical appearance of Hearlab.

2.0 Dimensions

2.1 Hearlab consists of the Stimulus Controller (SC), the Electrode Processor (EP), and a personal computer (PC).

2.2 The SC has the following dimensions: H: 12 cm, W: 32 cm, D: 29 cm.

2.3 The SC weights 2.7 kg.

2.4 The EP has the following dimensions: H: 5.5cm, W: 23 cm, D: 10.5 cm.

2.5 The EP weights 0.5 kg.

2.6 The PC dimensions and weight are not specified and subject to change.

3.0 Color Picture

3.1 Picture of the SC:



3.2 Picture of the EP:

