

INTRODUCING REAL-EAR MEASUREMENTS FONIX 8000 HEARING AID TEST SYSTEM



Software version 2.12 includes:

- Single or dual probe microphones with an advanced integrated earhook design
- Visible Speech
- Improved Digital Speech signal for testing open fit hearing aids
- New MOD NAL fitting formula for mild to moderate losses (based on NAL-NL1)

Compatible with the FONIX NOAH Real-ear Module!

Order the Real-ear Option with the 8120 Polar Plot sound chamber for the most advanced hearing aid testing on the market!

